

# MOUSE HYBRIDOMA PRODUCING IGG3 ANTIBODIES AGAINST HUMAN TIGIT (WUCAM)

---

[Boles, Kent, Colonna, Marco](#)

[Poranki, Deepika](#)

T-008862

Mouse IgG3 was generated by immunizing BALB/c mice with WUCAM<sup>FLAG</sup>-Baf3 cells. This can be used to study T cell-dependent antibody responses.

**Publication:** [A novel molecular interaction for the adhesion of follicular CD4 T cells to follicular DC](#)